


# CAPP *Harmony*

## SEROLOGICAL PIPETTES

**CAPP *Harmony* Serological Pipettes** are made from transparent medical-grade polystyrene sterilized by gamma irradiation. Each lot is tested for endotoxins (C.5 Eu/ML) (ISO T0993-11) and Cytotoxicity (ISO 10993-5) and are guaranteed non-haemolytic (ISO 10993-4). Ideal for use within tissue culture applications.

### Features:

- Designed for accurate and rapid fluid dispensing, these polystyrene pipettes are sterile, pyrogen-free and cotton-plugged for a variety of liquid handling applications
- The sharp, black, highly legible permanent graduations prevent ambiguity and guarantee precise reading and accurate dispensing to +1-2%
- Pipettes are supplied bulk (25 pipettes per bag) or individually wrapped in fibre-free colour-coded paper and plastic combination that is easy to open. See the ordering information below
- Each pipette comes with an effective filter barrier to prevent overflow of liquid into pipette
- Reverse graduations show the volume removed or remaining on all pipettes larger than 1mL.
- Generous negative graduations provide additional pipetting volume.
- The uniform mouth piece fits all major pipettors and is colour-coded for ease of identification



## Ordering Information

Cat. No.	Description	Volume	Package	Colour
SP-1-C	CAPPHarmony Serological pipette	1mL	individually wrapped, 500 pcs.	●
SPB-1-C	CAPPHarmony Serological pipette	1mL	bulk, 25 pcs./bag, 40 bags/case, (1000 pcs.)	●
SP-2-C	CAPPHarmony Serological pipette	2mL	individually wrapped, 500 pcs.	●
SPB-2-C	CAPPHarmony Serological pipette	2mL	bulk, 25 pcs./bag, 40 bags/case, (1000 pcs.)	●
SP-5-C	CAPPHarmony Serological pipette	5mL	individually wrapped, 200 pcs.	●
SPB-5-C	CAPPHarmony Serological pipette	5mL	bulk, 25 pcs./bag, 20 bags/case (500 pcs.)	●
SP-10-C	CAPPHarmony Serological pipette	10mL	individually wrapped, 200 pcs.	●
SPB-10-C	CAPPHarmony Serological pipette	10mL	bulk, 25 pcs./bag, 20 bags/case (500 pcs.)	●
SP-25-C	CAPPHarmony Serological pipette	25mL	individually wrapped, 150 pcs.	●
SPB-25-C	CAPPHarmony Serological pipette	25mL	bulk, 25 pcs./bag, 8 bags/case (200 pcs.)	●
SP-50-C	CAPPHarmony Serological pipette	50mL	individually wrapped, 100 pcs.	●
SPB-50-C	CAPPHarmony Serological pipette	50mL	bulk, 25 pcs./bag, 4 bags/case (100 pcs.)	●
SP-100-C	CAPPHarmony Serological pipette	100mL	individually wrapped , 50 pcs.	●